

Education
The Science of Hearing Aid Fitting

How Speechmapping Can Facilitate and Guide Solution-Centric Treatment

Presenter — David J. Smriga, M.A.
Senior Audiology Consultant
Audioscan Education

Kentucky Academy of Audiology
August 6, 2016

www.audioscan.com USA 800-265-2282 319-288-2218 info@audioscan.com
Copyright 2016 Electronic Design Incorporated. All rights reserved.

Any of these been on your mind lately?

Office of Science and Technology Policy

The National Academies of
SCIENCES • ENGINEERING • MEDICINE

OTC/DTC "Hearing Aids" ←

Audioscan® Education
Copyright 2016 Electronic Design Incorporated. All rights reserved.

Feelin' The Squeeze?

Delivery Options

Direct-to-Consumer

Lower Margins

Pricing Pressure

Perceived Value of Professional Care

Au.D.

Audioscan® Education
Copyright 2016 Electronic Design Incorporated. All rights reserved.

Key Part of NAS 12 Recommendations

- Recommendation 2:
 - Develop and promote measures to assess and improve quality of hearing health care services
 - Align and promote **best practices** and core competencies across the continuum of hearing health care, and implement mechanisms to insure widespread adherence
 - Research, develop and implement a set of quality metrics and measures to evaluate hearing health care services with the end goal of improving hearing and **communication focused patient outcomes**

Audioscan® Education
Copyright 2016 Pyramond Design Incorporated. All rights reserved.



Best Practices Guidelines – AAA 2006

(Modified by AuDNet Patient Care Excellence Program)

- **Provider Guarantee**
 - Auditory Assessment
 - Auditory Needs Assessment
 - Non-Auditory Needs Assessment
 - Hearing Instrument Selection
- * {
- Fitting & Verification
 - Hearing Instrument Orientation
 - Counseling & Follow-Up AR
 - Outcome Assessment
 - **Patient Care Audit**


*** Communication Focused Patient Outcomes**

Audioscan® Education
Copyright 2016 Pyramond Design Incorporated. All rights reserved.





Normal Hearing Hearing Loss




Adults with early-stage age-related hearing loss (right) show decreased activation of the hearing portion of the brain compared with normal hearing age-matched adults (left).

<http://www.hearingreview.com/2015/05/researchers-discover-brain-reorganizes-hearing-loss/>

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Normal Hearing Hearing Loss



Adults with mild age-related hearing loss (right) show brain reorganization in which hearing portions of their brain are recruited for processing visual patterns. This is not seen in age-matched adults with normal hearing (left).

Credit: Anu Sharma

<http://www.hearingreview.com/2015/05/researchers-discover-brain-reorganizes-hearing-loss/>

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Question. . .

- Can a rewired brain be rewired again?

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

INNOVATIVE APPROACHS

The brain's habit of rewiring itself takes time, an important concern in good fittings.

HABRAT: Hearing Aid Brain Rewiring Accommodation Time

STUART GATEHOUSE, PHD, AND MEAD C. KILLION, PHD

Every practicing dispenser knows that it may take a considerable amount of time for the listener to "get used to" new hearing aids. But do fundamental and beneficial changes in the wiring of the brain itself help the auditory rehabilitation process? Data presented here indicate it is sometimes impossible to properly evaluate the benefit of a particular hearing aid fitting for weeks or even months, which in turn implies that extraordinary judgment is required when the "wear it awhile and you'll get used to it" vs. "let me recast the hearing aid" decision must be made by the dispenser.

© PSC HR

S.C.: So I'll be happy with my hearing aids in another 3-6 weeks?
Dr. Abonzo: Maybe, maybe not. Before going through that much uncertainty, it would be a good thing to have your dispenser do a "real ear" measurement to make absolutely certain your hearing aids are giving you the right amount of high frequency gain. Perhaps they need to be readjusted.

S.C.: So after my hearing aids have been readjusted, I'll be happy in 3-6 more weeks and not toss them into my drawer drawer?
Dr. Abonzo: Unless the problem is that the hearing aids are painfully loud

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Hearing Instruments Volume 44, Number 10, 1992

Weeks	SHAPED RESPONSE (% Correct)	FLAT RESPONSE (% Correct)
0	75	75
2	76	74
4	75	73
6	78	75
8	79	75
10	81	75
12	82	75

5. Gatehouse S. The time course and magnitude of peripheral acclimatization to frequency responses: Evidence from monaural fitting of hearing aids. J Acoust Soc Am 92:1256-1268, 1992.

Audioscan® Educa
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Michael Merzenich on plasticity:

- Plasticity exists from cradle to grave
- Radical improvements in cognitive function are possible even in the elderly
- **Practicing a new skill under the right conditions** can change millions if not billions of connections between nerve cells in our brain maps

Doidge, N., "The Brain That Changes Itself" Penguin Books, 2007

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Kevin Munro's work

studying plasticity and acclimatization

- The ABR amplitude has been found to be greater in ears that had been aided, than in ears with the same degree of hearing impairment that had not been aided.
- The brain will only reorganize if speech is amplified to new, higher levels than the person previously experienced.

Audi: THE HEARING JOURNAL SEPTEMBER 2010 • VOL. 63 • NO. 9

Copyright 2016 Pyramind Design Incorporated. All rights reserved.



Two Key Take-aways

- The brain will only reorganize if speech is amplified to new, higher levels than the person previously experienced. (Munro)
- Practicing a new skill under the right conditions is required to change our brain maps. (Merzenich)

Audioscan® Education

Copyright 2016 Pyramind Design Incorporated. All rights reserved.



Aural Rehabilitation (AR)

“... the reduction of hearing-loss-induced deficits of function, activity, participation, and quality of life through sensory management, instruction, perceptual training and counseling.”

Boothroyd, A., "Aural Rehabilitation: What is it and Does it Work?"

Trends in Amplification, 2007 Jun; 11(2): 63-71

Audioscan® Education

Copyright 2016 Pyramind Design Incorporated. All rights reserved.




Treatment
It's a matter of perspective




Today's Techno-centric Model

What? COMMONLY USED MODEL




Jessen, D., Altidis, P., "ADA and IAA: Aural Rehabilitation in Private Practice"
ADA webinar recording at <http://audiologist.org/ada-and-iaa-aural-rehabilitation-in-private-practice>
Audiocan® Education
Copyright 2016. All rights reserved.




Proposed Solution-centric Model

What? PROPOSED MODEL



Jessen, D., Altidis, P., "ADA and IAA: Aural Rehabilitation in Private Practice"
ADA webinar recording at <http://audiologist.org/ada-and-iaa-aural-rehabilitation-in-private-practice>
Audiocan® Education
Copyright 2016. All rights reserved.




How Speechmapping Can Facilitate and Guide Solution-centric Treatment



Horizontal lines for notes

The Basic Speechmap Screen - Verifit




Audioscan Education Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Horizontal lines for notes


The Basic Speechmap Screen - Verifit



Audioscan Education Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Horizontal lines for notes



Killion, M., Mueller, H.G., "Twenty Years Later: The New Count-The-Dots Method", Hearing Journal, January 2010, pp 10-15



Child's Name: _____
 Test Date: _____
 NOTES: _____

INSTRUCTIONS

Step 1:

- Calculate the child's 5-frequency PTA on the 200, 400, 800 Hz.
- For scores 1-16 dB PTA, 4 dB from 200 to 200 Hz, 10 dB from 400 Hz.
- For scores 17-30 dB PTA, fit an oblique line as close as possible, connecting the reference frequencies of the fitting.

Step 2:

- Place the added SII value associated with the RTAB measured with a 20 dB SPL speech input on panel A relative to the corresponding PTA on the X-axis.
- Place the added SII value associated with the RTAB measured with a 40 dB speech input on panel B relative to the corresponding PTA on the X-axis.

INTERPRETATION:

- The solid line on each panel indicates the average added SII values for children aged birth to 4 years.
- The dashed lines indicate the mean confidence intervals around the mean.
- If the child's added SII values fall within the shaded areas, then the child's hearing aid fitting is a statistically acceptable fit for her PTA.


Audioscan® Education
 Copyright 2014 PedAMP Design Inc.


The Institute of Modern Otology Products Publishing of Hearing Aid Fitting Protocol Version 1.0, Revised 1
 ©2014 The Institute of Modern Otology Products Publishing of Hearing Aid Fitting


From "PedAMP" @ www.dslio.com




A Bit About Rules









DSL - Adult




NAL-NL2



NAL-NL2



NAL-NL2 REAR Results Using 5 Different Manufacturers' Programming Software



Sanders, J., Stoody, T., Weber, J., Meuller, H., "Manufacturers' NAL-NL2 Fittings Fail Real Ear Verification" Hearing Review, March 2015 24-32

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Pre-Fitting in the Test Box



RECD Applications - Reminder

- Convert audiometry (e.g., threshold and UCL) from HL to SPL near the eardrum
- Convert test box measurements of hearing aid output to estimated real-ear aided response (S-REM)

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



RECD Key Points

- For infants and children, RECD values vary from child to child and over time
- Measured RECD values can be different than average RECD values, even if average values are age specific – including adults



WRECD

Wideband Real-ear-to-coupler difference


- Difference in dB [across **200 - 12500 Hz**] between SPL measured in real-ear and in a **0.4cc** coupler, produced by transducer generating same wideband input signal.



WRECD

2cc coupler issue: Signals beyond 8KHz get buried in microphone noise floor


- ~ 14 dB greater output in 0.4cc coupler than 2cc coupler due to smaller volume
- Brings output above the measurement mic noise floor in HF
- System automatically converts 0.4cc RECD to HA-1 RECD when required






How do we measure RECD/WRECD ?

1. Calibrate the RECD transducer

WRECD: Verifit 2



RECD: Verifit VF-1, Axiom, RM500SL




Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



You are ready to measure the real ear response!

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Ideal End Result – VF2



WRECD

Max: 99.97% 120

RECD coupling: F100 to


Presentation: Single tone
Format: Graph
Scale (dB): SPL
Age: ALBC


Calibration: WRECD: 120dB
Last: Wed Apr 22 6:57 am
Interval: Weekly
Expires: Mon Apr 27 12:00 am


Prescribe real ear

Legend
Manufacturer calibration
Real ear
Basic WRECD
Avg WRECD

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.








Binaural test box



Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

New couplers and TRIC adapter


- New putty-less coupling for RIC instruments
- Slide-on couplers for increased speed.



Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Binaural test box

- Run simultaneous tests to compare programming or to test for binaural features.
- Instruments now positioned as worn on ear.



As-worn positioning of instruments

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.

Binaural Speechmap®

- Symmetrical, controlled acoustic space allows us to ensure same stimulus applied to each hearing instrument
- Use to test pairs, or to compare two different instruments!




Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Binaural Speechmap®

- Choose L+R link button in Speechmap dual view screen
- Tests will run simultaneously using one set of test parameters.



Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.




Calibrated Speech Signal Fitting Options

- Speech Standard 1 (the "Carrot Passage")
- Speech Standard 2
- ISTS

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.




A Speechmapping/Auditory Training Strategy



The Notion of “Goal” Fitting & “Starter” Fitting

- Definition of “Goal” Fitting:
 - Brings speech sounds back into the listening range for as broad a range of frequencies as possible:
 - Delivers an acceptable SII result
- Definition of “Starter” Fitting:
 - Settings that the patient feels comfortable enough with to start their new listening experience
- Difference between the two:
 - Practice!


Audiocan® Education
Copyright 2016. All rights reserved.



Calibrated Speech Signal Fitting Options

- Speech Standard 1 (the “Carrot Passage”)
- Speech Standard 2
- ISTS

Audiocan® Education
Copyright 2016. All rights reserved.




Initial Aided Audibility Verification

Yes or No – Is This A “Goal” Fitting?

- Set instrument to the most experienced user setting with the fitting software
- Select the fitting formula you wish to use
- Program the hearing instrument (First Fit)
- Place pre-programmed hearing instrument in the test box:
 - Session is re-loaded with RECD if available
- Run Speechmap at 65dB input in TEST 2 and obtain long-term average
 - Did you hit the target? **Yes or No**
 - Do you have an acceptable aided SII? **Yes or No**



Example of an Initial Aided Audibility Result



Speechmap® fitting protocol

- Test 1:
 - Set input to 50dB STD speech
 - This will show you a new target for that input level
 - Goal: To adjust the overall gain of the aid so that the middle line of the aided speech banana hits the target values indicated




Speechmap Fitting Protocol


- Test 3:
 - Input: MPO Sweep
 - Goal: To adjust the MPO of the aid so that the blue dots hit the output target, or are about 3dB below the UCL asterisks



Fitting Test Three: MPO



Final Fitting Summary




Some Brain Training Tools

Computer Based	Apps	Booklet	Group
•LACE* (neurotone.com) •Read My Quips* (sensesynergy.com) •Seeing & Hearing Speech # (sensimetrix.com)	•Hear Coach*## (itunes.apple.com) •Music Based Mobile Auditory Training Game*## (orca-us.info)	Cut To The Chase*## (cuttothechasehearingpros.com)	•G.R.O.U.P. AR* (idainstitute.com) •ACE – Active Communication Enhancement* https://shrs.uq.edu.au

* = Adult
= School age children

Audioscan® Education
Copyright 2016 Pyramont Design Incorporated. All rights reserved.



VERIFICATION CONSIDERATIONS


OPEN FIT TECHNOLOGIES


 Education
The Science of Hearing Aid Fitting


What Is An “Open” Instrument

- “Open” should NOT be defined by hearing aid design
 - Thin Tube
 - RIC
- “Open” defined by the way ANY hearing aid is coupled to the ear:
 - If coupling does not occlude the ear, it is an “Open” fitting

Audioscan® Education
Copyright 2016 Pyramont Design Incorporated. All rights reserved.



Minimal Occlusion



Lybarger S. Earmolds. In: Katz J, ed. Handbook of Clinical Audiology, 3rd edition. Baltimore: Williams and Wilkins, 1985: 885-910.

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Method For Determining "Open" Fit

- Record 65dB Speech (Test 1) with only probe tube in ear (REUR)
- Repeat 65dB Speech (Test 2) with aid now present, but muted (REOR)
- If two results are essentially the same, it is an "Open" fitting
- If Test 2 shows occlusion (i.e., loss of ear canal resonance), it is a "BTE" fitting

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



If an Open Fitting Is Present

- Pre-fit is off the table
 - There is no such thing as an "Open" coupler
- All fitting measures need to be done on the ear
- "Concurrent" calibration (equalization) should not be used


Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Impact of sound calibration method

For open fitting:

- Outflow from ear canal received by reference mic.
 - Lowers speaker output
- Resulting measured output will be reduced



Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Open Fitting Recommendation

- Use 'stored' equalization instead of 'concurrent' equalization
 - Disables reference mic.
 - Outflow from ear canal will not affect speaker level and your measurement
- If you cannot disable reference mic, move it away from ear canal



Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Effect of equalization approach




Figure 11. REAR findings for an OC fitting using concurrent equalization (reference microphone on) versus stored equalization (reference mic off)

Mueller, H.G., Ricketts, T.A., "Open Canal Fittings" Ten Take Home Tips", Hearing Journal, v9:11, Nov. 2006

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



An Open-Fit Speechmap Verification Protocol

1. Select "Open" in Instrument menu
2. Place probe tube / hearing aid on ear per typical
3. Turn OFF hearing aid.
4. Click on test signal and store 'equalization' when prompted
5. Turn ON hearing aid
6. Run Test 1 at 50dB calibrated speech to adjust gain
7. Run Test 2 at 65dB calibrated speech to adjust compression
8. Run Test 3 to set MPO

Optional:

1. Turn aid off and run TEST 4 at 65dB calibrated speech
2. Display Test 2 and Test 4 results as a counseling tool







FIGURE 5: The pink shaded area is the ear drum SPL "speech banana" for 65dB speech input measured at the probe tip with the open-fit hearing aid turned OFF. The green shaded area is the ear drum SPL "speech banana" with the same hearing aid turned ON. The difference between the two indicates where amplification has reached the eardrum.





Multiple choice




Audiocan Education
The Science of Hearing Aid Fitting

What is the core value proposition associated with an audibility-based counseling strategy?

- A. To restore normal hearing
- B. To make speech sound comfortable
- C. To bring speech sounds back into your listening range
- D. To make things sound a little bit louder than they do now


Audiocan® Education
Copyright 2016. All rights reserved.



Define an open fitting

- A. A fitting where the physical presence of the hearing instrument changes the acoustic properties of the ear canal
- B. A fitting where the physical presence of the hearing instrument does not change the acoustic properties of the ear canal
- C. A thin-tube or RIC hearing aid design
- D. A hearing aid or ear mold with a vent

Audiocan® Education
Copyright 2016. All rights reserved.



When should a .4cc coupler be used?

- A. To test RIC products
- B. To measure the gain of ITE hearing aids
- C. To measure aided output
- D. To effectively measure output for frequencies beyond 8KHz

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Case study discussion

Education
The Science of Hearing Aid Fitting

Case #1:

- *New patient*
- *Has never worn hearing aids before*
- *Hearing loss is mild-mod. S/N and symmetrical*
- *No thresholds greater than 70dBHL*

QUESTION: What is your initial counseling going to direct the patient to expect from your treatment?

QUESTION: How are you going to explain this concept to your patient?

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Case #2:

- *New patient*
- *Has never worn hearing aids before*
- *Hearing loss is mod.-severe S/N and symmetrical*
- *Thresholds of 80-85dB in the 3-8KHz range*

QUESTION: Are you going to direct your patient to expect anything different from your treatment than the previous patient?

Audioscan® Education
Copyright 2016 Plymouth Design Incorporated. All rights reserved.



Education
The Science of Hearing Aid Fitting

Thank You!

info@audioscan.com

audioscan
Hearing Instrument Fitting Systems

www.audioscan.com USA 800-265-2288 319-288-3318 info@audioscan.com
Copyright 2016 Plymouth Design Incorporated. All rights reserved.
